



## Stirrers

### iStir UNO+ - Slim High Capacity Stirrer

- iStir UNO +** Slim (29 mm only) single position Stirrer with 10 ltr capacity  
Brushless DC motor for maintenance free long life  
Set speeds from 200 RPM to 2200 RPM and Timer upto 999 min  
Pulse mode for CW & CCW direction change for thorough mixing  
Microprocessor controlled with gradual acceleration  
Corrosion resistant Stainless-Steel Top surface  
IP 44 class ingress protection (spill safe)  
Last run Memory feature  
Standard Accessories included:  
Supplied with 1 High Strength Stir bar (35 mm)  
Optional Accessories  
**B4** Stir Bar: 35 mm length



### iStir MS10 - Magnetic Stirrer

- iStir MS10** Noiseless and powerful motor - can stir upto 10 Lt media  
Digital display for ease of operation  
Speed setting from 200 to 2200 RPM  
Maximum capacity of 10 litres with large base for vessel  
Timer setting upto 999 minutes & infinite mode  
Programmable Pulse mode feature for vigorous mixing  
IP 42 protected design  
Stainless Steel top surface for corrosion resistance  
Closed design allows easy cleaning & avoids penetration of liquids  
**Standard Accessories included:**  
1 high strength Magnetic Stir Bar (25 x 8 mm)  
**Optional Accessories**  
**MS01-B1** Stir Bars (25 x 8 mm): 100 units per Pack



[CLICK HERE FOR MORE LAB EQUIPMENTS & SUPPLIES](#)